



Product Datasheet

EIF2C2 antibody (orb12795)


Description	Rabbit polyclonal antibody to EIF2C2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic peptide of human Argonaute-2(NP_036286.2).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	97kDa
Uniprot ID	Q9UKV8
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse brain tissue.



Immunohistochemical staining of mouse liver tissue.



Immunohistochemical staining of human tissue.