



Product Datasheet

Catenin-beta (phospho-Tyr489) antibody (orb127950)


Description	Rabbit Polyclonal to Catenin-beta (phospho-Tyr489).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P35222
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from C...



Western blot analysis of extracts from c...



Immunohistochemical analysis of formalin...