



Product Datasheet

Dynamin-1 (phospho-Ser774) antibody (orb127982)

Description	Rabbit Polyclonal to Dynamin-1 (phospho-Ser774).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q05193
Dilution Range	IHC-P: 1:50-1:100, IF/ICC: 1:100-1:500, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from m...



Immunohistochemical analysis of formalin...