



Product Datasheet

CRMP-2 (phospho-Ser522) antibody (orb127986)

Description	Rabbit Polyclonal to CRMP-2 (phospho-Ser522).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q16555
Dilution Range	WB: 1:500-1:3000, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from 3...



Western blot analysis of extracts from H...