



Product Datasheet

CIRBP antibody (orb12800)

Description	Rabbit polyclonal antibody to CIRBP
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human CIRBP(NP_001271.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	18kDa
Uniprot ID	Q14011
Dilution Range	WB: 1:500-1:1000, IP: 1:20-1:50, RIP: 1:20-1:50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U2OS cell...



Western blot analysis of extracts of var...