


Product Datasheet

PEK/PERK (phospho-Thr981) antibody (orb128012)

Description	Rabbit Polyclonal to PEK/PERK (phospho-Thr981).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9NZJ5
Dilution Range	IHC-P: 1:50-1:100, ELISA: 1:40000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...