



Product Datasheet

GSK3B antibody (orb12803)


Description	Rabbit polyclonal antibody to GSK3B
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthetic peptide of human GSK3 β (NP_001139628.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	46kDa
Uniprot ID	P49841
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:100, IF/ICC: 1:100-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Western blot analysis of extracts from n...



Western blot analysis of extracts of var...