



Product Datasheet

FGFR2 antibody (orb128062)


Description	Rabbit Polyclonal to FGFR2.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P21802
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from A...



Western blot analysis of extracts from H...



Immunohistochemical analysis of formalin...