



Product Datasheet

Gastrin Receptor antibody (orb128099)


Description	Rabbit polyclonal antibody to Gastrin Receptor
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P01350
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, IF/ICC: 1:100-1:500, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from N...



Western blot analysis of extracts from K...



Immunohistochemical analysis of formalin...