


Product Datasheet

GRIN2B (phospho-Ser1303) antibody (orb128142)

Description	Rabbit polyclonal antibody to GRIN2B phospho-Ser1303)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q13224
Dilution Range	WB: 1:500-1:3000, ELISA: 1:40000
Expiration Date	12 months from date of receipt.



Western
blot
analysis of
extracts
from M...