

Product Datasheet

IL22RA1 Antibody (orb1281567)

Description	Rabbit polyclonal antibody to IL22RA1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Fusion protein of human IL22RA1
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
Uniprot ID	Q8N6P7
Dilution Range	IHC: 1/25-1/100; ELISA: 1/2000-1/5000
Expiration Date	12 months from date of receipt.