

Product Datasheet

COL25A1 Antibody (orb1282420)

Description	Rabbit polyclonal antibody to COL25A1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Synthetic peptide of human COL25A1
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
MW	65 kDa
Uniprot ID	Q9BXS0
NCBI	NP_942014
Dilution Range	WB: 1/200-1/1000; IHC: 1/30-1/150; ELISA: 1/5000-1/10000
Expiration Date	12 months from date of receipt.