

Product Datasheet

C3/C3b Antibody (orb1282511)

Description	Rabbit polyclonal antibody to C3/C3b
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Synthetic peptide of human C3/C3b
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
MW	70 kDa
Uniprot ID	P01024
NCBI	NP_000055
Dilution Range	WB: 1/200-1/1000; IHC: 1/25-1/100; ELISA: 1/5000-1/10000
Expiration Date	12 months from date of receipt.