



Product Datasheet

Keratin 15 antibody (orb128263)


Description	Rabbit Polyclonal to Keratin 15.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	N-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P19012
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from C...



Western blot analysis of extracts from c...



Immunohistochemical analysis of formalin...