



Product Datasheet

Keratin 8 antibody (orb128271)

Description	Rabbit Polyclonal to Keratin 8.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P05787
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from S...



Immunohistochemical analysis of formalin...