

Product Datasheet

ARHGDIB Antibody (orb1282797)

Description	Rabbit polyclonal antibody to ARHGDIB
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Synthetic peptide of human ARHGDIB
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
MW	23 kDa
Uniprot ID	P52566
NCBI	NP_001166
Dilution Range	WB: 1/500-1/2000; IHC: 1/25-1/100; ELISA: 1/5000-1/10000
Expiration Date	12 months from date of receipt.