

Product Datasheet

RBFOX3 Antibody (orb1282901)

Description	Rabbit polyclonal antibody to RBFOX3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthetic peptide of human RBFOX3
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
Uniprot ID	A6NFN3
NCBI	NP_001076044
Dilution Range	IHC: 1/25-1/100; ELISA: 1/5000-1/10000
Expiration Date	12 months from date of receipt.