


Product Datasheet

MAK (phospho-Tyr159) antibody (orb128306)

Description	Rabbit Polyclonal to MAK (phospho-Tyr159).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P20794
Dilution Range	WB: 1:500-1:3000, ELISA: 1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from H...