

Product Datasheet

XKR3 Antibody (orb1283068)

Description	Rabbit polyclonal antibody to XKR3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthetic peptide of human XKR3
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
MW	53 kDa
Uniprot ID	Q5GH77
NCBI	NP_787074
Dilution Range	WB: 1/500-1/2000; ELISA: 1/2000-1/5000
Expiration Date	12 months from date of receipt.