

Product Datasheet

SLC5A11 Antibody (orb1283306)

Description	Rabbit polyclonal antibody to SLC5A11
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthetic peptide of human SLC5A11
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
Uniprot ID	Q8WWX8
NCBI	NP_443176
Dilution Range	IHC: 1/25-1/100; ELISA: 1/2000-1/5000
Expiration Date	12 months from date of receipt.