

Product Datasheet

DCAF6 Antibody (orb1283507)

Description	Rabbit polyclonal antibody to DCAF6
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthetic peptide of human DCAF6
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
Uniprot ID	Q58WW2
NCBI	NP_001017977
Dilution Range	IHC: 1/50-1/200; ELISA: 1/2000-1/5000
Expiration Date	12 months from date of receipt.