

Product Datasheet

IL5RA Antibody (orb1283658)

Description	Rabbit polyclonal antibody to IL5RA
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Synthetic peptide of human IL5RA
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
MW	48 kDa
Uniprot ID	Q01344
NCBI	NP_000555
Dilution Range	WB: 1/200-1/1000; IHC: 1/50-1/200; ELISA: 1/1000-1/2000
Expiration Date	12 months from date of receipt.