

## Product Datasheet

### IL31 Antibody (orb1283659)

|                            |                                                                                          |
|----------------------------|------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to IL31                                                       |
| <b>Species/Host</b>        | Rabbit                                                                                   |
| <b>Reactivity</b>          | Human                                                                                    |
| <b>Conjugation</b>         | Unconjugated                                                                             |
| <b>Tested Applications</b> | ELISA, IHC                                                                               |
| <b>Immunogen</b>           | Synthetic peptide of human IL31                                                          |
| <b>Preservatives</b>       | pH 7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol                                        |
| <b>Storage</b>             | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| <b>Note</b>                | For research use only                                                                    |
| <b>Isotype</b>             | Immunogen-specific rabbit IgG                                                            |
| <b>Clonality</b>           | Polyclonal                                                                               |
| <b>Uniprot ID</b>          | <a href="#">Q6EBC2</a>                                                                   |
| <b>NCBI</b>                | <a href="#">NP_001014358</a>                                                             |
| <b>Dilution Range</b>      | IHC: 1/50-1/200; ELISA: 1/2000-1/5000                                                    |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                          |