


Product Datasheet

APOBEC3B antibody (orb12845)

Description	Rabbit polyclonal antibody to APOBEC3B
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human APOBEC3B (NP_004891.4).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	46kDa
Uniprot ID	Q9UH17
Entrez	9582
Dilution Range	WB: 1:200-1:1000, IHC-P: 1:20-1:200, IP: 1:20-1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...