



Product Datasheet

SOCS3 antibody (orb128450)


Description	Rabbit polyclonal antibody to SOCS3.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P06748
Dilution Range	WB: 1:500-1:3000, IF/ICC: 1:100-1:500, IHC-P: 1:50-1:100, ELISA: 1:60000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from H...



Western blot analysis of extracts from J...



Immunohistochemical analysis of formalin...