

Product Datasheet

SLC2A11 Antibody (orb1284625)

Description	Rabbit polyclonal antibody to SLC2A11
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Synthetic peptide of human SLC2A11
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
Uniprot ID	Q9BYW1
NCBI	NP_110434
Dilution Range	IHC: 1/100-1/300; ELISA: 1/2000-1/10000
Expiration Date	12 months from date of receipt.