

Product Datasheet

GFRA2 Antibody (orb1284636)

Description	Rabbit polyclonal antibody to GFRA2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Synthetic peptide of human GFRA2
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
MW	52 kDa
Uniprot ID	O00451
NCBI	NP_001486
Dilution Range	WB: 1/500-1/2000; IHC: 1/10-1/50; ELISA: 1/1000-1/2000
Expiration Date	12 months from date of receipt.