



Product Datasheet

Trk B (phospho-Tyr515) antibody (orb128471)


Description	Rabbit polyclonal antibody to Trk B and is also known as GP145-TrkB whose expression can be detected...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q16620
Dilution Range	WB: 1:500-3000, ELISA: 1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from N...



Western blot analysis of extracts from K...



Immunohistochemical analysis of formalin...