

Product Datasheet

SMU1 antibody (orb12851)

Description	Rabbit polyclonal antibody to SMU1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human SMU1(NP_060695.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50kDa
Uniprot ID	Q2TAY7
Dilution Range	WB: 1:200-1:2000, IHC-P: 1:20-1:200, IF/ICC: 1:10-1:100
Expiration Date	12 months from date of receipt.