



## Product Datasheet

### PLD2 (phospho-Tyr169) antibody (orb128539)

|                            |                                                     |
|----------------------------|-----------------------------------------------------|
| <b>Description</b>         | Rabbit Polyclonal to PLD2 (phospho-Tyr169).         |
| <b>Species/Host</b>        | Rabbit                                              |
| <b>Reactivity</b>          | Human, Mouse, Rat                                   |
| <b>Conjugation</b>         | Unconjugated                                        |
| <b>Tested Applications</b> | IHC, WB                                             |
| <b>Immunogen</b>           | Internal                                            |
| <b>Concentration</b>       | 1mg/ml                                              |
| <b>Note</b>                | For research use only                               |
| <b>Clonality</b>           | Polyclonal                                          |
| <b>Uniprot ID</b>          | <b>O14939</b>                                       |
| <b>Dilution Range</b>      | WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:10000 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                     |



Western blot analysis of extracts from J...



Immunohistochemical analysis of formalin...