



Product Datasheet

60S Ribosomal Protein L10 antibody (orb128646)

Description	Rabbit Polyclonal to 60S Ribosomal Protein L10.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P27635
Dilution Range	WB: 1:500-1:3000, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from K...



Western blot analysis of extracts from K...