


Product Datasheet

SLC5A5 antibody (orb12865)

Description	Rabbit polyclonal antibody to SLC5A5
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein of human SLC5A5 (NP_000444.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	68kDa
Uniprot ID	Q92911
Entrez	6528
Dilution Range	WB: 1:500-1:2000, IF/ICC: 1:20-1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of mou...