

Product Datasheet

CCNE1 (Phospho-T395) antibody (orb12866)

Description	Rabbit polyclonal antibody to CCNE1 (Phospho-T395)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthetic phosphorylated peptide around T395 of human Cyclin E1(NP_001229.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	55kDa
Uniprot ID	P24864
Dilution Range	WB: 1:500-1:2000, IF/ICC: 1:50-1:100
Expiration Date	12 months from date of receipt.