

Product Datasheet

SERTAD3 Antibody (orb1286846)

Description	Rabbit polyclonal antibody to SERTAD3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Fusion protein of human SERTAD3
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
Uniprot ID	Q9UJW9
Dilution Range	IHC: 1/200-1/400; ELISA: 1/5000-1/10000
Expiration Date	12 months from date of receipt.