

Product Datasheet

ITGB3BP Antibody (orb1288157)

Description	Rabbit polyclonal antibody to ITGB3BP
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Fusion protein of human ITGB3BP
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
Uniprot ID	Q13352
Dilution Range	IHC: 1/50-1/200; ELISA: 1/1000-1/5000
Expiration Date	12 months from date of receipt.