

Product Datasheet

AIF1 Antibody (orb1288183)

Description	Rabbit polyclonal antibody to AIF1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Full length fusion protein
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
MW	17 kDa
Uniprot ID	P55008
Dilution Range	WB: 1/200-1/1000; IHC: 1/50-1/200; ELISA: 1/1000-1/2000
Expiration Date	12 months from date of receipt.