



Product Datasheet

Tyrosine Hydroxylase (phospho-Ser19) antibody (orb128835)

Description	Rabbit polyclonal antibody to Tyrosine Hydroxylase (phospho-Ser19).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	N-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P07101
Dilution Range	IHC-P: 1:50-1:100, IF/ICC: 1:100-1:500, ELISA: 1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from 2...



Immunohistochemistry analysis of formali...