



Product Datasheet

Thyroid Hormone Receptor alpha antibody (orb128839)

Description	Rabbit Polyclonal to Thyroid Hormone Receptor alpha.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	N-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P10827
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from SKOV3 cells.



Immunohistochemical analysis of formalin-fixed paraffin-embedded (FFPE) tissue sections.