


Product Datasheet

TR-beta1 (phospho-Ser142) antibody (orb128840)

Description	Rabbit Polyclonal to TR-beta1 (phospho-Ser142).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P10828
Dilution Range	WB: 1:500-1:3000, ELISA: 1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from 2...