



Product Datasheet

TIMP1 antibody (orb128843)

Description	Rabbit polyclonal antibody to TIMP1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P01033
Dilution Range	WB: 1:500-1:3000, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from 2...



Western blot analysis of extracts from c...