



Product Datasheet

TNNI3K antibody (orb128856)

Description	Rabbit Polyclonal to TNNI3K.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q59H18
Dilution Range	IHC-P: 1:50-1:100, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from c...



Immunohistochemical analysis of formalin...