

Product Datasheet

c-Myc Tag Antibody (orb1288616)

Description	Rabbit polyclonal antibody to c-Myc Tag
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthetic peptide: C-EQKLISEEDL conjugated with KLH.
Preservatives	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Immunogen-specific rabbit IgG
Clonality	Polyclonal
MW	55 kDa
Uniprot ID	P01106
Dilution Range	WB: 1/10000-1/50000; ELISA: 1/5000-1/10000
Expiration Date	12 months from date of receipt.