

Product Datasheet

M6PR Antibody (orb1288672)

Description	Mouse monoclonal antibody to M6PR
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IHC
Immunogen	Purified recombinant fragment of human M6PR (AA: 27-277) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	31kDa
Entrez	4074
Dilution Range	IHC: 1/200 - 1/1000; ICC: 1/50 - 1/200; RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.