



Product Datasheet

TSC22D1 antibody (orb128874)


Description	Rabbit Polyclonal to TSC22D1.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	N-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q15714
Dilution Range	WB: 1:500-1:3000, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from m...



Immunohistochemical analysis of formalin...



Confocal immunofluorescence analysis of ...