



Product Datasheet

Tubulin alpha antibody (orb128878)


Description	Rabbit Polyclonal to Tubulin alpha.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P68363, Q71U36
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, IF/ICC: 1:100-1:500, ELISA: 1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from M...



COLO and HepG2 cells using Tubulin alpha...



Western blot analysis of extracts from 3...