

Product Datasheet

SARS-CoV-2-ACE2 Antibody (orb1289005)

Description	Mouse monoclonal antibody to SARS-CoV-2-ACE2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Purified recombinant fragment of SARS-CoV-2-ACE2 expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Dilution Range	ELISA: 1/10000
Expiration Date	12 months from date of receipt.