



Product Datasheet

CCR10 antibody (orb1289030)


Description	Mouse monoclonal antibody to CCR10
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	Purified recombinant fragment of human CCR10 (AA: extra mix) expressed in E. Coli.
Target	CCR10
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Form/Appearance	Purified antibody in PBS with 0.05% sodium azide
Concentration	1.0 mg/ml
Storage	4°C; -20°C for long term storage
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	38.4kDa
Entrez	2826
Dilution Range	WB 1/500 - 1/2000, FCM 1/200 - 1/400, ELISA 1/10000
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
HL-60 cells ...



Flow
cytometric
analysis of
MOLT4 cells
...



Flow
cytometric
analysis of
THP-1 cells ...