

Product Datasheet

YB1 (phospho-Ser102) antibody (orb128919)

Description	Rabbit Polyclonal to YB1 (phospho-Ser102).	Image not found or type unknown
Species/Host	Rabbit	Western blot
Reactivity	Human, Mouse, Rat	analysis of
Conjugation	Unconjugated	extracts
Tested Applications	WB	from H...
Immunogen	Internal	
Concentration	1mg/ml	
Note	For research use only	
Clonality	Polyclonal	
Uniprot ID	P67809	
Dilution Range	WB: 1:500-1:3000, ELISA: 1:40000	
Expiration Date	12 months from date of receipt.	