


Product Datasheet

Ferrous Iron Microplate Assay Kit (orb1289377)

Description	Ferrous Iron Microplate Assay Kit
Reactivity	All
Range	5 µmol/L - 500 µmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Colorimetric
Sample Types	Serum, Plasma, Urine, Saliva, Tissue extracts, Cell lysate, Other biological fluids
Expiration Date	Please enquire.



Representative
Standard
Curve