

Product Datasheet

Mouse Immunoglobulin Isotyping ELISA Kit (orb1289380)

Description	Mouse Immunoglobulin Isotyping ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Colorimetric
Sample Types	Cell Culture Supernates, Ascites, Purified Antibodies Preparations
Expiration Date	Please enquire.